


# MEDDP-519

19" Slim Medical LCD Monitor



CE FC



## Features

- 19" LCD Panel with LED Backlight
- Compact and Terse Design
- VGA x 1, DVI x1
- 5 wire resistive touch screen
- CE/ FCC Class B Passed 18, EN60601-1:2006
- Low Power Consumption

## Specifications

### Display

Size	19"
Resolution	1280 x 1024
Max. Colors	16.7M
Brightness (Typ.)	350 nits
Viewing Angle (CR=10)(Typ.)	170°(H)/160°(V)
Contrast Ratio (Typ.)	1000:1
Touch Screen	5 wire resistive touch screen (optional)

### System

OS Support	Windows® XP, Windows® 7, Windows® 10
DICOM Support	Optional

### I/O

Audio	Audio-in x 1
Serial Port	USB x 1 for touch controller
Signal	VGA x 1, DVI x 1
Speaker	1W speaker x 1 on back side

### Mechanical and Environmental

Power Consumption	Full Loading :35.5Watts
Power Requirement	DC12V / 5A power input
Operating Temperature	5°C ~ 40°C(41°F ~ 104°F)
Storage Temperature	-10°C ~ 60°C(14°F ~ 140°F)
Storage Humidity	10%~90%@35°C, non-condensing
Gross Weight	5.8kg (12.79lb)
Net Weight	4.9kg (10.8lb)
Dimension	421.1 x 350.5 x 48.1 mm
Package Dimension (W x D x H)	490 x 500 x 155 mm
Certifications	CE: EN 60601-1-2:2007(V3.0) , EN 60601-1:2006/ A1:2013/A12:2014(V3.1) FCC: Part 18 Class B, Part 15 Class B IEC 60601-1: 2012 (V3.1)

## Ordering Information

- » **MEDDP-519ANN-A1-1010**  
19" LCD Monitor, 350 nits,VGAx1,DVI x1, DC 12V Adapter
- » **MEDDP-519ARN-A1-1010**  
19" LCD Monitor, 350 nits, Resistive Touch, VGAx1,DVI x1,DC 12V Adapter
- » **MEDDP-519TND-A1-1010** (DICOM)  
19" LCD Monitor, 350 nits, VGAx1.HDMIx1.DPx1, DC 12V Adapter, DICOM Function
- » **MEDDP-519TRD-A1-1010** (DICOM)  
19" LCD Monitor, 350 nits, Resistive Touch VGAx1.HDMIx1.DPx1, DC 12V Adapter, DICOM Function